

HPLC LABORATORY SYSTEMS


separpur PCS

PEPTIDE PURIFICATION SYSTEM

separpur PCS chromatographs present compact universal and reliable units for laboratory scale prep HPLC. Modification of basic system is possible in broad range to tailor the final product according customer needs.

separpur PCS chromatographs are designed to operate with high pressure, high performance columns having inner diameter in a range 10 mm - 100 mm. They consist of:

- Pulse less pump **separatrix** PP03
- Photometric (UV or UV-VIS) detector
- Prep 6 port loop injector
- Carousel fraction collectors **separflow** FFC 96
- Keyboard and display
- Clarity chromatostation or/and with ECOMAC MS software
- Option accessories according demand



Type	Pump	Detector	Software	Additional
SEPARPUR PCS 400S	SEPARTRIX PP03SA	Flash 10 DAD 600	Clarity Lite Chromatostation	
SEPARPUR PCS 400SG				Set 3 gradient valves
SEPARPUR PCS 400SGC				Fraction collector